



AMBULATORY AND IN-HOSPITAL EEG

Improved Patient Care


- Continuous data acquisition
- Patient can move freely
- Continuous impedance checking
- EEG storage up to 2000Hz per channel

IP Camera Switching

- Switch cameras on the fly as patients move
- Change recording cameras from acquisition station, review station, or Sentinel Monitoring station

Sentinel Monitoring

- Monitor multiple patients from one computer
- Get immediate visual and auditory notifications for patient events and seizure detections
- Control camera functions
- See EEG traces and trends



Arc Apollo[™] EEG

AT-HOME EEG SOLUTION

Easy to Use

- Up to 96 hours of EEG recording on a single battery charge
- Programmable patient event buttons
- Microphone for voice events
- Rechargeable, nonproprietary lithium ion batteries

Q-Video Mobile 3

- HD video with auto-switching infrared

Designed for Patient Comfort

- Compact, 1 lb 3 oz. (539 g) and rugged
- Backpack, harness, and head mount options

